ITALIANO

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ENGLISH

Instructions for use of prosthetic parts for Modus implant and compatible

Abstract: This "COMPATIBLE" prosthetic line is designed to be used for the preparation of dental implants (internal hexagon and external hexagon type) of the Modus system, produced by Erreceiffe.

Prosthetic parts Class IIb:  0425

HEALING ABUTMENT (Class IIb - made of titanium gr.5)

Remove the cover screw of the implant. Clean the head of the implant with a saline solution, ensuring that it is free from any presence of bone or gum tissue. Measure the height of the soft tissue and then choose the height of the healing abutment to be used. Screw the healing abutment. Recommended permanence duration 15-20 days. Tightening Torque of the healing abutment: 15Ncm

TITANIUM ABUTMENT for MUA - implants - EWB (Class IIb - made of titanium gr.5)

It is equipped with the fixing screw with which we connect it to the implant. It is suitable for screwed prostheses (bridges and single crowns). Due to its non-rotating hexagonal system and its conical shape it is used to support both single crowns and bridges. It is possible to shorten and modify it, and then model it directly on the implant. The EWB type is suitable for the realization of bars thanks to the threaded portion. Tightening Torque of the locking screw: 32Ncm

TITANIUM ABUTMENT (type STANDARD o T-BASE) (Class IIb - made of titanium gr.5)

It consists of two parts and is mainly used for total fixed prostheses, for bars or when aesthetic is not the only requisite for the success of the case. It is recommended when we focus on proper maintenance of oral hygiene. It allows you to perform a fixed screwed prosthesis even in the presence of marked divergences. Tightening Torque of the locking screw: 32Ncm

ANGLED ABUTMENT (Class IIb - made of titanium gr.5)

It consists of a single piece. It is used when the implants present significant non-parallelism, which is difficult to correct with straight abutments. Tightening Torque of the locking screw: 32Ncm

TITANIUM ABUTMENT (type MUA Abument) (Class IIb - made of titanium gr.5)

It consists of a single piece. It is mainly used in total fixed prosthesis. It is particularly suitable for cases with a limited intra-maxillar space, with excellent aesthetics possibilities. Tightening Torque of the locking screw: 32Ncm

BALL ATTACHMENT (Class IIb - made of titanium gr.5)

Suitable for anchoring overdentures prostheses in cases where the divergence of the implants does not exceed 10°. For use with the retentive cap made of nylon, usually located in the dedicated steel case (available on the market). Screwing torque: 32Ncm

 **GOLD ABUTMENT (GOLD type) (Class IIb - the castable portion is made of POM, the screw is made of titanium gr.5, the base is made of alloy of gold PLATINOR N).** The concept of this abutment is that of the castable. The advantage is that the connection to the implant is already pre-formed and does not have to be created by fusion. The base made of gold alloy is available with two types of connection to the implant: with hexagon (non-rotating) and without hexagon (rotating). The base is connected to the castable portion through a cylindrical connection with a diameter of 3.5 mm. The castable portion is fully adjustable by the technician. The version with the hexagon is ideal for single tooth screwed or for the construction of abutments. The version without hexagon is suitable for screwed bridges, for substructures connected directly to the implants or for overdentures bars. Tightening torque of the locking screw: 32Ncm

MULTI UNIT ABUTMENT (MUA) STRAIGHT AND ANGLED (Class IIb - Made in Titanium Alloy - TI-6AL-4V)

The angled abutments are equipped with a fixing screw and with a screwing-rod. The Multi-unit Abument system is designed to provide the maximum versatility and to improve the aesthetic result of multiple screwed prostheses. The low vertical height of the straight abutments and the angled abutments of 17 and 30 degrees allow an appropriate prosthesis even in the most difficult situations. They also offer the possibility of organizing the implants with a maximum angulations of 45°. The reference clinical solution is the "All-on-4", developed to maximize the use of the available bone portion and to employ the technology of the immediate loading. Using only four implants in edentulous jaws, this solution takes advantage of the inclination of the posterior implants to provide a secure and optimal prosthetic support for a bridge, even with minimum bone volume. Tightening torque of the locking screw of straight MUA: === ⚡ 32Ncm
Tightening torque of the locking screw angled MUA of 17 ° and 30 °: === ⚡ 15Ncm

Prosthetic parts Class I: 

CASTABLE ABUTMENT (type UCLA) (Class I - made of meltable plastic POM. The fixing screw is made of titanium grade 5)

It allows the maximum personalization of the prosthetic product. Two versions are available: with hexagon (non-rotating) and without hexagon (rotating). The type with the hexagon is indicated for individual teeth screwed directly on the implant or for the construction of abutments. The type without hexagon is suitable for screwed bridges and substructures connected to the implants. Also available without hex for MUA abutments. Tightening torque of the locking screw: 32Ncm

SCAN ABUTMENT (Class I - made of PEEK, The fixing screw is made of titanium grade 5)

Used to detect, with the help of scanners, the exact DIGITAL position of the implant. It is used as a mean of communication between the laboratory and the dentist to build abutments through CAD software. Tightening torque of the locking screw: 32Ncm

ANALOG (Class I - made of titanium gr.5)

Connect the analog of the implant to the transfer with the fixing screw. In this phase it is necessary to place considerable attention to the fixing maneuvering, in order to ensure the precise coupling between the implant analogue and the transfer. To avoid torsions to the transfer incorporated in the impression, it is recommended to lock down the analog during this phase and not to over-tighten the screw.

TRANSFER (Class I - made of titanium gr.5)

Place the precision transfer on the implant with its locking screw which, in this case, is longer than the transfer. This allows the screw to come out of the impression case and to be unscrewed easily. Detect the impression exclusively with an individual perforated tray, in correspondence of the implant, to allow the identification of the fixing screw. Wait for the perfect setting of the impression material, unscrew the screw with the special driver and remove the impression. In this case the transfer remains embedded in the impression material without losing its position.

EASY WAY BAR (Class I - Made of Titanium gr.5)

Kit for 4-5-6 implants. The blocks are screwed onto the threaded abutments and act as a hook for the telescopic arms. Type 1 blocks are used for the frontal area, while type 2 blocks are used in the posterior area. The arms are indicated to connect to type 1 or 2 blocks by means of a spherical connection. They can have different lengths, be straight or arched / male or female and they are used as a connection between the implants / abutments.

CAUTION:

- The medical devices are distributed NON-STERILE. If provided, sterilize before their use for 20 minutes at 134 ° C
- Exclusive dental use only
- Always verify the integrity of the packaging

WARNING: DO NOT REUSE a device classified as DISPOSABLE. Even if is not evident, it may be mechanically deformed (thread or hexagon) or it may be contaminated.

More detailed information on the use of the medical device are given in the specific guidelines. If you do not have it, ask to your distributor or download it directly from our website in the section "DOCUMENTS".

	FABBRICANTE-MANUFACTURER ERRECIFFE srl Via V. Emanuele II, 68 - 24036 PONTE SAN PIETRO (BG) - ITALIA	www.errecliffe.com E-mail: info@errecliffe.com
---	---	--

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